

MC045-1

Mesa de exploración Ginecológica 2 niveles



Design

Simple structure, practical; Whole stainless steel material meets the needs of the surgical room. This bed is designed for gynecological diagnosis, operation and parturition.



Material

Baosteel from Fortune Global 500. Panasonic Robotic ensure 360° full smooth welding.



Leather

Flame retardant, yellow resistance, mildew resistance, bending resistance, 100,000 times wear resistance.



Sponge

Flame retardant, 60mm thickness, high-density foam, density 28kg/m3.



Lithotomy Rod

Lithotomy Rod, adjustable in height and different angles.



Hand Bar

Granular surface to increase friction and prevent hand slip.



Stainless Steel Basin

Removable stainless steel basin, convenient to take out.



Feet Cover

4 pcs Anti-slip bed feet covers.



Parámetros Técnicos

Longitud Externa	1820mm
Ancho Externo	590mm
Altura	800mm

Configuración Técnica

✓ Varilla de litotomía	1 set	✓ Palanca de mando	2sets
✓ Barras para manos	1 set	✓ Bandeja de acero inoxidable	1ud
✓ Tapizado antideslizante	4uds		

