

MCB2A12

Tabla de Inmovilización espinal



Material

This product is blow molded by high density plastic polyethylene, strong and durable, not easy to aging.




Function

Compatible with X-ray, MRI, CT transparent, convenient for examination of patient on board, reduce the possible hurt to unload the patient from the board.



Handle

The handle whole is evenly opened around, which can be used for many people to lift, carry and lift at the same time.



Float

Hard structure, easy to continue CPR and cardiac compression rescue in the process of transport, large displacement, the overall volume of 0.04m3, can float an adult in normal temperature water, greatly reducing the difficulty of water rescue personnel, convenient for faster rescue of the injured.



Bandages

This product is equipped with three bandages as standard. The optional kit includes neck bracket, head immobilizer and spine fix plate.



Control de calidad importante

- * Disponible con 10 pines
- * Radiografía translúcida
- * Compatible con la mayoría de los inmovilizadores de cabeza
- * De una pieza
- * Impermeable a todos los fluidos corporales y fácil de descontaminar.
- * Es el mejor tipo flotante
- * Color: naranja / amarillo
- * Material: polietileno de alta densidad (HDPE)
- * USFDA

Parámetros técnicos

Tamaño completo L*A*A

1840*450*55mm

